## DIAQUICK hCG Combo Cassette

*A rapid test for the qualitative detection of human chorionic gonadotropin (hCG) human urine, serum and plasma samples*

### REF

#### Content

**Z01405CE** - 30 tests individually packed in foil pouches, disposable pipette (30x REF Z01405B)  
 - 1 package insert

**For professional in vitro diagnostic use only**

### INTENDED USE

The DIAQUICK hCG Combo Cassette is a rapid chromatographic immunoassay for the qualitative detection of human chorionic gonadotropin (hCG) in urine, serum or plasma, to aid in the early detection of pregnancy.

### DIAGNOSTIC IMPLICATIONS

Human chorionic gonadotropin (hCG) is a glycoprotein hormone secreted by the developing placenta shortly after fertilization. In normal pregnancy, hCG can be detected in serum, plasma and urine as early as 7 to 10 days after conception.<sup>1,2,3,4</sup> hCG levels continue to rise very rapidly, frequently exceeding 100 mIU/mL by the first missed menstrual period,<sup>2,3,4</sup> and peaking in the 100,000-200,000 mIU/mL range about 10-12 weeks into pregnancy. The appearance of hCG soon after conception, and its subsequent rapid rise in concentration during early gestational growth, make it an excellent marker for the early detection of pregnancy.

The DIAQUICK hCG Combo Cassette is a rapid test to qualitatively detect the presence of hCG in urine, serum and plasma specimens at the sensitivity of 25 mIU/mL. The test utilizes a combination of monoclonal and polyclonal antibodies to detect hCG. The test shows no cross-reactivity with the structurally related glycoprotein hormones hFSH, hLH and hTSH at physiological levels.

### TEST PRINCIPLE

The test uses two lines to indicate results. The test line utilizes a combination of antibodies including a monoclonal hCG antibody to detect hCG. The control line is composed of goat polyclonal antibodies and colloidal gold particles. The assay is conducted by adding specimen to the specimen well of the cassette and observing the formation of coloured lines. The specimen migrates via capillary action along the membrane to react with the coloured conjugate. Positive specimens react with the specific antibody-hCG-coloured conjugate to form a coloured line at the test line region of the membrane. Absence of this coloured line suggests a negative result. To serve as a procedural control, a coloured line will always appear in the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

### STORAGE

Store as packaged in the sealed pouch at room temperature or refrigerated (2-30 °C). The kit is stable until the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

### MATERIALS PROVIDED

- test cassettes
- disposable pipettes
- package insert

### MATERIALS REQUIRED BUT NOT PROVIDED

- specimen collection container
- timer

### REAGENTS

The test cassette contains anti-hCG particles and anti-hCG coated on the membrane.

### SAMPLE COLLECTION AND PREPARATION

#### Urine Assay

A urine specimen must be collected in a clean and dry container. A first morning urine specimen is preferred since it generally contains the highest concentration of hCG; however, urine specimens collected at any time of the day may be used. Urine specimens exhibiting visible precipitates should be centrifuged, filtered, or allowed to settle to obtain a clear specimen for testing.

#### Serum or Plasma Assay

Blood should be collected aseptically into a clean tube without anticoagulants to obtain serum or with anticoagulants to obtain plasma. K<sub>2</sub>EDTA, Na-citrate, K-citrate, Na-heparin, Li-heparin, Na-oxalate can be used as anticoagulants. Separate the serum or plasma from blood as soon as possible to avoid haemolysis. Use clear non-haemolysed specimens when possible.

#### Specimen Storage

Specimens may be stored at 2-8 °C for up to 48 hours prior to testing. For prolonged storage, specimens may be frozen and stored below -20 °C. Frozen specimens should be thawed and mixed before testing.


### PRECAUTIONS

- For professional in vitro diagnostic use only. Do not use after the expiration date.
- The test should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test should be discarded according to local regulations.

### ASSAY PROCEDURE

**Allow the test cassette, specimen and/or controls to equilibrate to room temperature (15-30 °C) prior to testing.**

1. Bring the pouch to room temperature before opening it. Remove the test cassette from the sealed pouch and use it as soon as possible.
2. Place the test cassette on a clean and level surface. Hold the dropper vertically and **transfer 3 full drops of urine, serum or plasma** (approx. 120 µL) to the specimen well (S) of the test cassette, and then start the timer. Avoid trapping air bubbles in the specimen well (S). See the illustration below.
3. Wait for the coloured line(s) to appear. **Read the result at 3 minutes when testing a urine specimen, or at 5 minutes when testing a serum or plasma specimen.** Do not interpret results after 10 minutes. It is important that the background is clear before the result is read.

### INTERPRETATION OF RESULTS

**POSITIVE:** Two distinct coloured lines appear. One line should be in the control region (C) and another line should be in the test region (T). One line may be lighter than the other; they do not have to match. This result means a pregnancy is probable.

**NEGATIVE:** One coloured line appears in the control region (C). No apparent red or pink line appears in the test region (T). This means that a pregnancy is not likely.

**INVALID:** Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

### LIMITATIONS

1. The DIAQUICK hCG Combo Cassette is a qualitative test, therefore, neither the quantitative value nor the rate of increase in hCG can be determined by this test.
2. Very diluted urine specimens, as indicated by a low specific gravity, may not contain representative levels of hCG. If pregnancy is still suspected, a first morning urine specimen should be collected 48 hours later and tested.
3. Very low levels of hCG (less than 50 mIU/mL) are present shortly after implantation. However, because a significant number of first trimester pregnancies terminate for natural reasons, a test result that is weakly positive should be confirmed by retesting with a first morning urine or serum specimen collected 48 hours later.
4. This test may produce false positive results. A number of conditions other than pregnancy, including trophoblastic disease and certain non-trophoblastic neoplasms, including testicular tumours, prostate cancer, breast cancer, and lung cancer, cause elevated levels of hCG. Therefore, the presence of hCG in urine or serum specimens should not be used to diagnose pregnancy unless these conditions have been ruled out.
5. This test may produce false negative results. False negative results may occur when the levels of hCG are below the sensitivity level of the test. When pregnancy is still suspected, a first morning urine or serum specimen should be collected 48 hours later and tested. In case pregnancy is suspected and the test continues to produce negative results, see a physician for further diagnosis.
6. As with any assay employing mouse antibodies, the possibility exists for interference by human anti-mouse antibodies (HAMA) in the specimen. Specimens from patients who have received preparations of monoclonal antibodies for diagnosis or therapy may contain HAMA. Such specimens may cause false positive or false negative results.
7. This test provides a presumptive diagnosis for pregnancy. A confirmed pregnancy diagnosis should only be made by a physician after all clinical and laboratory findings have been evaluated.
8. This test reliably detects intact hCG up to 1,000,000 mIU/mL.

### QUALITY CONTROL

A procedural control is included in the test. A red coloured line appearing in the control region (C) is the internal positive procedural control. It confirms sufficient specimen volume and correct procedural technique. A clear background is an internal negative procedural control. If a background colour appears in the result window and interferes with the ability to read the test result, the result may be invalid. It is recommended that a positive hCG control (containing 25-250 mIU/mL hCG) and a negative hCG control (containing "0" mIU/mL hCG) be evaluated to verify proper test performance when a new shipment of test devices are received.

### EXPECTED VALUES

Negative results are expected in healthy non-pregnant women and healthy men. Healthy pregnant women have hCG present in urine, serum or plasma specimens. The amount of hCG will vary greatly with gestational age and between individuals. The DIAQUICK hCG Combo Cassette has a sensitivity of 25 mIU/mL, and is capable of detecting pregnancy as early as 1 day after the first missed menses.

## PERFORMANCE CHARACTERISTICS

### Accuracy

A multi-centre clinical evaluation was conducted comparing the results obtained using the DIAQUICK hCG Combo Cassette and another commercially available urine/serum/plasma membrane hCG test. The urine study included 413 specimens and both assays identified 296 negative and 117 positive results. The serum/plasma study included 200 specimens and both assays identified 141 negative and 59 positive results. The results demonstrated a > 99.0 % overall accuracy of the DIAQUICK hCG Combo Cassette when compared to the other membrane hCG test.

#### hCG Reference Method (urine)

Method	Other hCG Rapid Test		Total Results
	Positive	Negative	
DIAQUICK hCG Combo Cassette	Results		
	Positive	117	0
	Negative	0	296
<b>Total Results</b>		117	296

Relative Sensitivity: 100.0 % (96.9 % - 100 %) \*

Relative Specificity: 100.0 % (98.8 % - 100 %) \*

Accuracy: 100.0 % (99.1 % - 100 %) \*

\* 97.5 % Confidence Interval

#### hCG Reference Method (serum, plasma)

Method	Other hCG Rapid Test		Total Results
	Positive	Negative	
DIAQUICK hCG Combo Cassette	Results		
	Positive	59	0
	Negative	0	141
<b>Total Results</b>		59	141

Relative Sensitivity: 100.0 % (93.9 % - 100 %) \*

Relative Specificity: 100.0 % (97.4 % - 100 %) \*

Accuracy: 100.0 % (98.2 % - 100 %) \*

\* 97.5 % Confidence Interval

### Sensitivity and Specificity

The DIAQUICK hCG Combo Cassette detects hCG at concentrations of 25 mIU/mL or greater. The test has been standardized to the WHO 4<sup>th</sup> International Standard (75/589). The addition of LH (300 mIU/mL), FSH (1,000 mIU/mL), and TSH (1,000 µIU/mL) to negative (0 mIU/mL hCG) and positive (25 mIU/mL hCG) specimens showed no cross-reactivity.

### Precision

#### Intra-Assay

Within-run precision has been determined by using 10 replicates of four specimens containing 25 mIU/mL, 100 mIU/mL, 250 mIU/mL and 0 mIU/mL of hCG. The negative and positive values were correctly identified 100 % of the time.

#### Inter-Assay

Between-run precision has been determined by using the same four specimens of 25 mIU/mL, 100 mIU/mL, 250 mIU/mL and 0 mIU/mL of hCG in 10 independent assays. Three different lots of the DIAQUICK hCG Combo Cassette have been tested. The specimens were correctly identified 100 % of the time.

### Interfering Substances

The following potentially interfering substances were added to hCG negative and positive specimens.

Acetaminophen	20 mg/dL	Caffeine	20 mg/dL
Acetylsalicylic Acid	20 mg/dL	Gentisic Acid	20 mg/dL
Ascorbic Acid	20 mg/dL	Glucose	2 g/dL
Atropine	20 mg/dL	Haemoglobin	1 mg/dL
Bilirubin (serum/plasma)	40 mg/dL	Bilirubin (urine)	2 mg/dL
Triglycerides (serum/plasma)	1200 mg/dL		

None of the substances at the concentration tested interfered in the assay.

## REFERENCES

1. Batzer FR. Hormonal evaluation of early pregnancy. *Fertil. Steril.* 1980; 34(1): 1-13
2. Catt KJ, ML Dufau, JL Vaitukaitis. Appearance of hCG in pregnancy plasma following the initiation of implantation of the blastocyst. *J. Clin. Endocrinol. Metab.* 1975; 40(3): 537-540
3. Braunstein GD, J Rasor, H. Danzer, D Adler, ME Wade Serum human chorionic gonadotropin levels throughout normal pregnancy. *Am. J. Obstet. Gynecol.* 1976; 126(6): 678-681
4. Lenton EA, LM Neal, R Sulaiman Plasma concentration of human chorionic gonadotropin from the time of implantation until the second week of pregnancy. *Fertil. Steril.* 1982; 37(6): 773-778
5. Steier JA, P Bergsjö, OL Myking Human chorionic gonadotropin in maternal plasma after induced abortion, spontaneous abortion and removed ectopic pregnancy. *Obstet. Gynecol.* 1984; 64(3): 391-394
6. Dawood MY, BB Saxena, R Landesman. Human chorionic gonadotropin and its subunits in hydatidiform mole and choriocarcinoma. *Obstet. Gynecol.* 1977; 50(2): 172-181
7. Braunstein GD, JL Vaitukaitis, PP Carbone, GT Ross Ectopic production of human chorionic gonadotropin by neoplasms. *Ann. Intern Med.* 1973; 78(1): 39-45

## USED SYMBOLS

Symbol	Description
	Content
	Dispose of the tests and packaging appropriately

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## BIBLIOGRAPHIE

1. Batzer FR. Hormonal evaluation of early pregnancy. *Fertil. Steril.* 1980; 34(1): 1-13
2. Catt KJ, ML Dufau, JL Vaitukaitis. Appearance of hCG in pregnancy plasma following the initiation of implantation of the blastocyte, *J. Clin. Endocrinol. Metab.* 1975; 40(3): 537-540
3. Braunstein GD, J Rasor, H. Danzer, D Adler, ME Wade Serum human chorionic gonadotropin levels throughout normal pregnancy, *Am. J. Obstet. Gynecol.* 1976; 126(6): 678-681
4. Lenton EA, LM Neal, R Sulaiman Plasma concentration of human chorionic gonadotropin from the time of implantation until the second week of pregnancy, *Fertil. Steril.* 1982; 37(6): 773-778
5. Steier JA, P Bergsjö, OL Myking Human chorionic gonadotropin in maternal plasma after induced abortion, spontaneous abortion and removed ectopic pregnancy, *Obstet. Gynecol.* 1984; 64(3): 391-394
6. Dawood MY, BB Saxena, R Landesman. Human chorionic gonadotropin and its subunits in hydatidiform mole and choriocarcinoma, *Obstet. Gynecol.* 1977; 50(2): 172-181
7. Braunstein GD, JL Vaitukaitis, PP Carbone, GT Ross Ectopic production of human chorionic gonadotropin by neoplasms, *Ann. Intern Med.* 1973; 78(1): 39-45

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